

Positive Control - T-Cell
Anti-Human T Lymphocyte
Mouse MoAb (Monoclonal Antibody)

Package Insert

PC-T, 0.5 mL
SPC-T, 50 µL

INTENDED USE

For use as a control that reacts with human T-lymphocytes.
Use undiluted.

For Research Use Only.

SUMMARY/EXPLANATION

This mouse monoclonal antibody reacts positively with human T-lymphocytes. The material is a mouse ascitic fluid diluted in heat inactivated fetal bovine sera.

Warning: Animal products and human blood products have the potential to transmit disease. Take appropriate precautions.

Caution: This product contains <0.1% sodium azide. Sodium azide reacts with lead and copper plumbing forming highly explosive metal azides.

Dispose of according to local regulations.

INSTRUMENTS/EQUIPMENT
(Not Included with reagent)

≤-50°C freezer

INSTRUCTIONS FOR USE

1. Thaw just prior to aliquoting.
2. Aliquot into microcentrifuge tubes in volume suitable for single use.
3. Freeze aliquots immediately at ≤-50°C.
4. Thaw only as needed for use.

REACTIVITY

Determine reactivity in your laboratory's test system.

STORAGE, HANDLING, SHIPPING

Positive control for T Cells is shipped frozen, on dry ice. Promptly store the frozen material at ≤-50°C.



GTi DIAGNOSTiCS®
Good science starts with people.®

Positive Control - T-Cell

- FOR RESEARCH USE ONLY
- STORE AT ≤-50°C

20925 Crossroads Circle, Suite 200
Waukesha, WI 53186-4054 USA
(262) 754-1000 or 1-800-233-1843

Cat. NO. PC-T, SPC-T
Rev. 12 October 2007

www.gtidiagnostics.com