CATALOG #3275 Catalogue # 3275

For *in vitro* diagnostic use

Contenu:
50 SPECIMEN COLLECTION KITS
Each kit contains:
One urethral/conjunctival swab
One tube containing 0.7 mL
GEN-PROBE transport medium
Store at Room Temperature

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Diagnostic *"in vitro"* uniquement

Use only with the GEN-PROBE® PACE® Systems and AMPLIFIED™ Chlamydia Trachomatis Assay

**SPECIMEN
COLLECTION
KITS FOR URETHRAL OR
CONJUNCTIVAL SPECIMENS**

Use only with the GEN-PROBE PACE® Systems and AMPLIFIED Chlamydia Trachomatis Assay

For *in vitro* diagnostic use
Store at Room Temperature

**SPECIMEN
COLLECTION
KITS FOR URETHRAL OR
CONJUNCTIVAL SPECIMENS
CATALOG #3275**

Specimen Collection and Handling for Male Urethral Specimens

Patient should not have urinated for at least 1 hour prior to sample collection.

Instructions for Chlamydia trachomatis (PACE and AMPLIFIED) and/or Neisseria gonorrhoeae (PACE only)

1. Insert swab from collection kit 2 to 4 cm into urethra.
2. Once inserted, gently rotate swab clockwise using sufficient pressure to ensure swab comes into contact with all urethral surfaces. Allow swab to remain inserted for 2 to 3 seconds.
3. Withdraw swab.
4. Prepare swab for transport.


Specimen Collection and Handling for Chlamydia trachomatis Conjunctival Specimens (PACE Only)

1. If pus or discharge is present, use a sterile Dacron swab (not provided) to clean the area. Do not scrape the conjunctiva while cleaning the eye(s).
2. If both eyes are affected, swab the least affected eye first.
3. Thoroughly swab the lower then the upper conjunctiva 2 to 3 times each with the GEN-PROBE urethral/conjunctival swab.
4. Prepare swab for transport.

Preparation for Transport

1. Fully insert one swab into the GEN-PROBE transport tube.
2. Snap off shaft at score line or cut shaft to fit tube. Use care to avoid splashing of contents.
3. Cap tube tightly.
4. For best results, transport tube to the laboratory at 2° to 25°C. Store at 2° to 25°C and test within 7 days of collection.
5. If not tested within 7 days:
 - a. PACE: Vortex specimen, remove swab, and store at -20° to -70°C.
 - b. AMPLIFIED: Store at -20° to -70°C.**Do not remove swab.**

Note: Specimens collected with this system cannot be used for culture. Only GEN-PROBE urethral/conjunctival swabs should be used for specimen collection.

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