

Quik-ID Class I Recording Sheet RS

LOT 122110-QID

EXP 18 AUG 12



2012-08-18

NAME _____

ID# _____

BDT _____

TECH _____

DATE _____

<NC> (A) _____

PRA _____

| WELL | ETHN | GTI# | HLA | | | | | | BAF | 2xNC | CTF | SAMP | DIFF | INTR |
|----------|------|---------|-------|-------|-----|-------|-------|------|-------|---------|---------|------|-------|------|
| | | | | | | | | (B) | (C) | (D) | (E) | (F) | | |
| A1 | H | 140DFW | A24 | B35 | B61 | (-,6) | Cw3 | Cw4 | 0.893 | | | | | |
| B1 | C | 141CBC | A2 | A29 | B44 | B60 | (4,6) | Cw3 | Cw16 | 1.045 | | | | |
| C1 | C | 117CBC | A11 | A26 | B38 | B60 | (4,6) | Cw3 | Cw12 | 0.861 | | | | |
| D1 | AA | 15LB | A23 | A80 | B45 | B81 | (-,6) | Cw16 | Cw18 | 1.005 | | | | |
| E1 | C | 158CBC | A32 | A68 | B57 | B62 | (4,6) | Cw3 | Cw6 | 0.922 | | | | |
| F1 | C | 181DFW | A1 | A32 | B44 | B64 | (4,6) | Cw6 | Cw8 | 0.998 | | | | |
| G1 | C | 24HC | A24 | A32 | B7 | B48 | (-,6) | Cw7 | Cw8 | 0.921 | | | | |
| H1 | C | 160CBC | A11 | A29 | B7 | B56 | (-,6) | Cw1 | Cw7 | 0.837 | | | | |
| A2 | C | 112CBC | A2 | A24 | B18 | B27 | (4,6) | Cw7 | Cw7 | 1.017 | | | | |
| B2 | C | 194DFW | A3 | A11 | B27 | B60 | (4,6) | Cw1 | Cw3 | 1.060 | | | | |
| C2 | As | YH106 | A11 | A31 | B39 | B61 | (-,6) | Cw3 | Cw7 | 1.000 | | | | |
| D2 | As | 19HC109 | A31 | A33 | B35 | B46 | (-,6) | Cw1 | Cw4 | 1.008 | | | | |
| E2 | C | 22HC114 | A24 | A26 | B7 | B13 | (4,6) | Cw6 | Cw7 | 1.031 | | | | |
| F2 | U | PHL262 | A3 | A32 | B55 | B64 | (-,6) | Cw3 | Cw8 | 0.898 | | | | |
| G2 | As | CAS321 | A2 | A34 | B46 | B75 | (-,6) | Cw1 | Cw4 | 0.958 | | | | |
| H2 | AA | WAR374 | A2 | A23 | B45 | B50 | (-,6) | Cw6 | | 1.002 | | | | |
| A3 | AA | NSB267 | A29 | A68 | B51 | B62 | (4,6) | Cw3 | Cw16 | 1.005 | | | | |
| B3 | AA | ACD412 | A2 | A29 | B13 | B62 | (4,6) | Cw3 | Cw6 | 0.885 | | | | |
| C3 | C | BOG477 | A2 | A25 | B8 | B27 | (4,6) | Cw1 | Cw7 | 0.749 | | | | |
| D3 | AA | 50LB136 | A2 | A6601 | B7 | B53 | (4,6) | Cw4 | Cw7 | 0.932 | | | | |
| E3 | C | JMN461 | A1 | A3 | B57 | B57 | (4,-) | Cw6 | Cw18 | 1.106 | | | | |
| F3 | C | LER780 | A31 | A32 | B8 | B39 | (-,6) | Cw7 | Cw12 | 0.748 | | | | |
| G3 | AA | JKI546 | A2 | A23 | B42 | B53 | (4,6) | Cw6 | Cw17 | 0.798 | | | | |
| H3 | AA | 17LB634 | A1 | A6602 | B52 | B58 | (4,-) | Cw7 | Cw12 | 0.982 | | | | |
| A4 | AA | 73LB206 | A3 | A74 | B8 | B45 | (-,6) | Cw7 | Cw16 | 0.971 | | | | |
| B4 | AA | 64LB364 | A25 | A33 | B55 | B58 | (4,6) | Cw3 | Cw6 | 0.963 | | | | |
| C4 | C | TY636 | A2 | A30 | B7 | B56 | (-,6) | Cw1 | Cw7 | 0.986 | | | | |
| D4 | AA | 56LB758 | A23 | A33 | B41 | B47 | (4,6) | Cw7 | Cw8 | 1.071 | | | | |
| E4 | AA | 52LB017 | A30 | A80 | B18 | B78 | (-,6) | Cw2 | Cw16 | 0.992 | | | | |
| F4 | AA | 96LB831 | A30 | A36 | B58 | B65 | (4,6) | Cw6 | Cw8 | 1.061 | | | | |
| G4 | As | TCH077 | A1 | A2 | B37 | B76 | (4,6) | Cw4 | Cw6 | 0.986 | | | | |
| H4 | As | AHN822 | A26 | A30 | B51 | B62 | (4,6) | Cw3 | Cw14 | 0.944 | | | | |
| A5 | As | HII402 | A11 | A24 | B59 | B67 | (4,6) | Cw1 | Cw7 | 0.908 | | | | |
| B5 | C | AK342 | A1 | A69 | B8 | B73 | (-,6) | Cw7 | Cw15 | 0.788 | | | | |
| C5 | As | THI458 | A24 | A26 | B60 | B77 | (4,6) | Cw7 | Cw8 | 0.695 | | | | |
| D5 | AA | ICH467 | A29 | A68 | B49 | B71 | (4,6) | Cw3 | Cw7 | 0.801 | | | | |
| E5 | AA | NT927 | A1 | A2 | B35 | B52 | (4,6) | Cw4 | Cw12 | 0.876 | | | | |
| F5 | AA | 49LB684 | A3 | A33 | B47 | B63 | (4,-) | Cw6 | Cw14 | 0.833 | | | | |
| G5 | AA | JA118 | A33 | A74 | B53 | B72 | (4,6) | Cw2 | Cw4 | 0.882 | | | | |
| H5 | U | MAD981 | A2 | A2 | B54 | B60 | (-,6) | Cw1 | Cw7 | 0.833 | | | | |
| A,B,C,D6 | | POOL | NC | | | | | | | | | | | |
| E6 | | POOL | PC | | | | | | | | | | | |
| F6 | | None | NAW | | | | | | | | | | | |
| | | | CALCS | | | | | | | C=2x(A) | D=Bx(C) | | F=E-D | |



Quik-ID Class I Recording Sheet RS

Wells
Ag

LOT 122110-QID

wells containing antigen

| AG | # | AG | # |
|-------|----|-----|---|
| A1 | 6 | B7 | 5 |
| A2 | 12 | B8 | 4 |
| A3 | 5 | B13 | 2 |
| A11 | 5 | B18 | 2 |
| A23 | 4 | B27 | 3 |
| A24 | 6 | B35 | 3 |
| A25 | 2 | B37 | 1 |
| A26 | 4 | B38 | 1 |
| A28 | 0 | B39 | 2 |
| A29 | 5 | B40 | 0 |
| A30 | 4 | B41 | 1 |
| A31 | 3 | B42 | 1 |
| A32 | 5 | B44 | 2 |
| A33 | 5 | B45 | 3 |
| A34 | 1 | B46 | 2 |
| A36 | 1 | B47 | 2 |
| A43 | 0 | B48 | 1 |
| A66 | 0 | B49 | 1 |
| A6601 | 1 | B50 | 1 |
| A6602 | 1 | B51 | 2 |
| A68 | 3 | B52 | 2 |
| A69 | 1 | B53 | 3 |
| A74 | 2 | B54 | 1 |
| A80 | 2 | B55 | 2 |
| | | B56 | 2 |
| Cw1 | 8 | B57 | 2 |
| Cw2 | 2 | B58 | 3 |
| Cw3 | 12 | B59 | 1 |
| Cw4 | 7 | B60 | 5 |
| Cw5 | 0 | B61 | 2 |
| Cw6 | 11 | B62 | 4 |
| Cw7 | 17 | B63 | 1 |
| Cw8 | 6 | B64 | 2 |
| Cw12 | 4 | B65 | 1 |
| Cw13 | 0 | B67 | 1 |
| Cw14 | 2 | B70 | 0 |
| Cw15 | 1 | B71 | 1 |
| Cw16 | 5 | B72 | 1 |
| Cw17 | 1 | B73 | 1 |
| Cw18 | 2 | B75 | 1 |
| | | B76 | 1 |
| | | B77 | 1 |
| Bw4,6 | 22 | B78 | 1 |
| Bw4,- | 3 | B81 | 1 |
| Bw-,6 | 15 | B82 | 0 |

| WELL | GTI# | HLA-A | | HLA-B | | HLA-C | |
|------|---------|----------|--------|----------|----------|-------------|----------|
| A1 | 140DFW | 24:02 | | 35:12 | 40:02 | 03:06 | 04:01 |
| B1 | 141CBC | 02:01 | 29:02 | 40:01 | 44:03 | 03:04/06/09 | 16:01 |
| C1 | 117CBC | 11:01 | 26:01 | 38:01 | 40:01 | 03:04 | 12:03 |
| D1 | 15LB | 23:01 | 80:01 | 45:01 | 81:01 | 16:01 | 18:01/02 |
| E1 | 158CBC | 32:01 | 68:01 | 15:01 | 57:01 | 03:03 | 06:02 |
| F1 | 181DFW | 01:01 | 32:01 | 14:01 | 44:03 | 06:02 | 08:02 |
| G1 | 24HC | 24:02 | 32:01 | 07:02/06 | 48:01 | 07:02 | 08:03 |
| H1 | 160CBC | 11:01 | 29:02 | 07:02 | 56:01 | 01:02 | 07:02 |
| A2 | 112CBC | 02:01 | 24:02 | 18:01 | 27:05 | 07:01 | 07:04 |
| B2 | 194DFW | 03:01 | 11:01 | 27:05 | 40:01 | 01:02 | 03:04 |
| C2 | YH106 | 11:01 | 31:01 | 39:01 | 40:02 | 03:04 | 07:02 |
| D2 | 19HC109 | 31:01 | 33:03 | 35:01 | 46:01 | 01:02 | 04:01 |
| E2 | 22HC114 | 24:03 | 26:01 | 07:02 | 13:02 | 06:02 | 07:02 |
| F2 | PHL262 | 03:01 | 32:01 | 14:01 | 55:01 | 03:03 | 08:02 |
| G2 | CAS321 | 02:07 | 34:01 | 15:21 | 46:01 | 01:02 | 04:03 |
| H2 | WAR374 | 02:05 | 23:01 | 45:01/02 | 50:01/02 | 06:02 | |
| A3 | NSB267 | 29:02 | 68:01 | 15:01 | 51:01 | 03:04 | 16:02 |
| B3 | ACD412 | 02:01 | 29:02 | 13:02 | 15:01 | 03:04 | 06:02 |
| C3 | BOG477 | 02:01 | 25:01 | 08:01 | 27:05 | 01:02 | 07:01 |
| D3 | 50LB136 | 02:02/05 | 66:01 | 07:02 | 53:01 | 04:01 | 07:02 |
| E3 | JMN461 | 01:01 | 03:01G | 57:01 | 57:02 | 06:02 | 18:01 |
| F3 | LER780 | 31:01 | 32:01 | 08:01 | 39:01 | 07:01 | 12:03 |
| G3 | JKI546 | 02:01 | 23:01 | 42:01 | 53:01 | 06:02 | 17:01G |
| H3 | 17LB634 | 01:01 | 66:02 | 52:01 | 58:01 | 07:18 | 12:02 |
| A4 | 73LB206 | 03:01 | 74:01 | 08:01 | 45:01 | 07:01 | 16:01 |
| B4 | 64LB364 | 25:01 | 33:03 | 55:01 | 58:01 | 03:03 | 06:02/07 |
| C4 | TY636 | 02:01 | 30:01 | 07:02 | 56:01 | 01:02 | 07:02 |
| D4 | 56LB758 | 23:01 | 33:03 | 41:01 | 47:01 | 07:01G | 08:02 |
| E4 | 52LB017 | 30:01 | 80:01 | 18:01 | 78:01 | 02:02 | 16:01 |
| F4 | 96LB831 | 30:02 | 36:01 | 14:02 | 58:02 | 06:02 | 08:02 |
| G4 | TCH077 | 01:01 | 02:03 | 15:12 | 37:01 | 04:03 | 06:02/03 |
| H4 | AHN822 | 26:02 | 30:01 | 15:01 | 51:01 | 03:03 | 14:02 |
| A5 | HI1402 | 11:01 | 24:02 | 59:01 | 67:01 | 01:02 | 07:02 |
| B5 | AK342 | 01:01 | 69:01 | 08:01 | 73:01 | 07:01 | 15:05 |
| C5 | THI458 | 24:02 | 26:01 | 15:13 | 40:01 | 07:02 | 08:01 |
| D5 | ICH467 | 29:02 | 68:02 | 15:10 | 49:01 | 03:04 | 07:01G |
| E5 | NT927 | 01:01 | 02:01 | 35:01 | 52:01 | 04:01 | 12:02 |
| F5 | 49LB684 | 03:01 | 33:03 | 15:16 | 47:01 | 06:02 | 14:02 |
| G5 | JA118 | 33:03 | 74:01 | 15:03 | 53:01 | 02:02 | 04:01 |
| H5 | MAD981 | 02:03 | 02:07 | 40:01 | 54:01 | 01:02 | 07:02 |

| SYMBOL | ENGLISH | DEUTSCH | FRANÇAIS | GREEK | ITALIANO | PORTUGUÊS | ESPAÑOL |
|--------|---|--|--|---|--|--|--|
| RS | Recording Sheet | Protokollbogen | Tableaux de résultats | Φύλλο Καταγραφής | Foglio di lavoro | Tabela de resultados | Tabla de Resultados |
| NAME | Patient Name | Name des Patienten | Nom du patient | Όνομα Ασθενούς | Nome paziente | Nome do paciente | Nombre de paciente |
| ID # | Identification Number | Identifikationsnummer | N° identification du patient | Αριθμός Ταυτότητας | Numero di identificazione | Número de identificação | N° de identificación |
| HLA | HLA Type | HLA Typ | Type de système HLA | Τύπος HLA | Tipizzazione HLA | Tipo HLA | Tipo HLA |
| BDT | Bleed Date | Entnahmedatum | Date de prélèvement | Ημερομηνία αιμοληψίας | Data prelievo | Data da perda de sangue | Fecha de Extracción |
| <NC> | Mean of Negative Controls (A) | Durchschnittswert der negativen Kontrollen (A) | Moyenne des contrôles négatifs (A) | Μέσος Όρος Τυμών Αρνητικού Ορού Ελέγχου (A) | Media dei controlli negativi (A) | Média dos Controlos Negativos (A) | Media Controles Negativos (A) |
| DATE | Date Read | Datum | Date | Ημερομηνία | Data | Data da leitura | Fecha |
| TECH | Technologist | Untersucher | Technologue | Τεχνολόγος | Tecnico | Técnico | Analista |
| PRA | Percent Reactive Antibody | Panel Reaktiver Antikörper % | Pourcentage d'anticorps réactifs | Ποσοστό Αντιδρώντων Αντισωμάτων | Percentuale reattività anticorpi | Percentagem de anticorpo reactivo | Porcentaje Reactividad Ab |
| GTI # | GTI # | HLA Typ | N° GTI | GTI# | GTI# | GTI# | GTI# |
| ETHN | Ethnicity of HLA Source | Ethnische Ursprung | Ethnie de la source HLA | Εθνογένεση της πηγής των HLA | Sorgente etnica degli HLA | Etnia da fonte de HLA | Grupo Racial del origen de HLA |
| WELLS | # of wells containing antigen | Frequenz der Antigene | Nombre de puits avec antigène | Αριθμός θετικών πηκτών που περιέχουν το αντίγονο | N° di pozzetti contenenti l'antigene | N° de poços com antígeno | Número pocillo conteniendo Ag |
| AG | Antigen | Antigen | Antigène | Αντίγονο | Antigene | Antígeno | Antígeno |
| BAF | Background Adjustment Factor (B) 2 times NC value (C) | Background Adjustment Faktor (B) 2 fois la valeur NC (C) | Facteur d'ajustement de fond (B) 2 fois la valeur NC (C) | Παράγοντας Ρύθμισης Υποβάθρου (B) Η τιμή του αρνητικού ορού ελέγχου επί 2 (C) | Fattore aggiustamento Colorazione di fondo (B) Valore dei controlli neg. X 2 (C) | Factor de Ajustamento do Background (B) 2 vezes o valor NC (C) | Factor Ajuste del Cero (B) 2 veces el valor NC (C) |
| CTF | Cutoff (D) | Cutoff (D) | Coupure (D) | Όριο θετικότητας (D) | Cutoff (D) | Cut-off (D) | Cut Off (D) |
| SAMP | Sample (E) | Probe (E) | Echantillon (E) | Δείγμα (E) | Campione (E) | Amostra (E) | Muestra (E) |
| DIFF | Difference (F) | Differenz (F) | Différence (F) | Διαφορά (F) | Differenza (F) | Diferença (F) | Diferença (F) |
| WELL | Well | Vertiefung | Puits | Βόθρο | Pozzetto | Poço | Pocillo |
| INTR | Interpretation | Auswertung | Interprétation | Ερμηνεία | Interpretazione | Interpretação | Interpretación |
| NC | Negative Control | Negative Kontrolle | Contrôle négatif (A) | Αρνητικός Όρος ελέγχου | Controllo negativo | Controlo negativo | Control negativo |
| PC | Positive Control | Positive Kontrolle | Contrôle positif | Θετικός Όρος ελέγχου | Controllo positivo | Controlo positivo | Control positivo |
| NAW | No Antigen Well | Veriefung ohne Antigen | Puits sans antigène | Βόθρο χωρίς αντίγονο | Pozzetto senza antigene | Poço sem antígeno | N° Pocillo Antígeno |
| CALCS | Calculations | Berechnungen | Calculs | Υπολογισμοί | Calcoli | Cálculos | Calculos |
| MoAb | Monoclonal Antibody | Monoklonaler Antikörper | Anticorps monoclonal | Μονοκλωνικό Αντίσωμα | Anticorpo monoclonale | Anticorpo monoclonal | Anticuerpo Monoclonal |
| DIL | Dilute Before Use | Vor Gebrauch verdünnen | Diluée avec de l'eau distillée avant utilisation | Αραιώστε πριν την χρήση | Diluire con acqua deionizzata prima dell'uso | Diluir | Dilución |

