

LIFECODES Class I ID
LM1
LOT 090111-LM1

BEAD	DONOR	HLA-A		HLA-B		HLA-C	
102	316	A*2601	A*3301	B*1402	B*4201	C*0802	C*1701
103	634	A*0101	A*6602	B*5201	B*5801	C*0718	C*1202
104	458	A*2402	A*2601	B*1513	B*4001	C*0702	C*0801
111	461	A*0101	A*0301	B*5701	B*5702	C*0602	C*1801
112	116	A*2301	A*8001	B*4501	B*8101	C*1601	C*1801
113	247	A*0101	A*6901	B*0801	B*3508	C*0701	C*1203
114	205	A*1101	A*2901	B*3505	B*4601	C*0102	C*0401
115	941	A*3002	A*6802	B*1516	B*4403	C*0701	C*1402
117	209	A*2402	A*3001	B*1302	B*4001	C*0602	C*0702
118	262	A*0301	A*3201	B*1401	B*5501	C*0303	C*0802
119	234	A*2601	A*6602	B*4501	B*5601	C*0102	C*0602
121	498	A*3301	A*3402	B*3501	B*5301	C*0401	
122	127	A*3002		B*1402	B*7801	C*0802	C*1601
125	106	A*1101	A*3101	B*3901	B*4002	C*0304	C*0702
126	531	A*3402	A*8001	B*1801	B*5301	C*0202	C*0602
127	267	A*2902	A*6801	B*1501	B*5101	C*0304	C*1602
128	792	A*2301	A*3601	B*0705	B*5801	C*0401	C*0701
129	012	A*2402	A*3101	B*3501	B*5201	C*0401	C*1202
130	517	A*1101		B*1502	B*5401	C*0102	C*0801
132	031	A*0201	A*3001	B*1503	B*5001	C*0210	C*0401
133	484	A*3001	A*3201	B*3508	B*7301	C*0401	C*1505
134	541	A*2501	A*3201	B*1801	B*4901	C*0701	C*1203
135	181	A*4301		B*1503		C*1802	
136	119	A*0301	A*3101	B*3501	B*5101	C*0401	C*1502
137	088	A*3001		B*4201	B*5001	C*0602	C*1701
138	397	A*0308	A*3402	B*1510	B*5703	C*0304	C*0701
139	307	A*0101	A*0201	B*0801	B*3906	C*0701	C*0702
141	311	A*0301	A*2902	B*5101	B*5601	C*0102	C*1502
142	372	A*0201	A*2301	B*3503	B*4006	C*0401	C*1502
143	398	A*0202	A*3002	B*1801	B*7801	C*0202	C*1601
144	626	A*3601	A*7401	B*5001	B*5802	C*0602	
145	412	A*0201	A*2902	B*1302	B*1501	C*0304	C*0602
146	325	A*0201	A*2601	B*2705	B*3801	C*0102	C*1203
147	490	A*2301	A*6602	B*0702	B*4410	C*0401	C*0702
148	402	A*1101	A*2402	B*5901	B*6701	C*0102	C*0702
149	544	A*2402	A*3303	B*1501	B*3701	C*0304	C*0602
150	114	A*2403	A*2601	B*0702	B*1302	C*0602	C*0702
151	876	A*6802	A*7401	B*0801	B*1503	C*0202	C*0701
152	467	A*2902	A*6802	B*1510	B*4901	C*0304	C*0701
154	730	A*0101	A*0201	B*0702	B*3701	C*0602	C*0702
155	812	A*2402	A*6801	B*4901	B*8201	C*0302	C*0701
156	758	A*2301	A*3303	B*4101	B*4701	C*0701	C*0802
157	117	A*2402	A*3201	B*3801	B*4801	C*0801	C*1203
161	479	A*2501	A*3303	B*4001	B*5301	C*0304	C*0401
162	770	A*0101	A*3303	B*0702	B*5703	C*0701	C*0702
163	477	A*0201	A*2501	B*0801	B*2705	C*0102	C*0701
164	295	A*3303		B*1516	B*5801	C*0302	C*1402
165	765	A*2402	A*6801	B*3901	B*5501	C*0303	C*0702
166	908	A*2501	A*6801	B*2702	B*2705	C*0102	C*0202
167	814	A*0101	A*6601	B*1402	B*4102	C*0802	C*1701

Note: In cases where several alleles are possible only the most common is represented in this chart. Refer to the worksheet supplied with each kit for details.