



LIFECODES Class II IDv2

LM2Q

LOT 072811-LM2QR



2012-07-25

EXP: 25 JUL 12

BEAD	ETHN	DONOR	HLA			
106	AA	814	DR11	DR1303	52	DQ7
113	C	373				DQ2
114	AA	126				DQ2
115	AA	031				DQ4
117	U	013				DQ5
118	AA	634				DQ6 DQ6
119	C	699				DQ7
121	U	839				DQ8
122	U	012				DQ9
124	As	109				DQ4 DQ9
125	C	461				DQ5 DQ9
126	AA	484	DR11	DR17	52	DQ2 DQ7
127	AA	467	DR8	DR15	51	DQ6 DQ7
128	As	106				DQ6 DQ8
129	AA	560				DQ4 DQ7
130	C	373	DR7		53	DQ2
132	AA	876	DR1	DR13	52	DQ5 DQ6
133	AA	490	DR9	DR18	52,53	DQ2 DQ4
135	C	590	DR7	DR17	52,53	DQ2 DQ2
136	AA	812	DR9		53	DQ2
137	C	328	DR1	DR4	53	DQ5 DQ7
138	AA	653	DR1	DR11	52	DQ5 DQ6
139	AA	941	DR1			DQ5
141	As	458	DR4	DR12	52,53	DQ7 DQ8
142	C	780	DR12	DR17	52	DQ2 DQ7
143	As	077	DR16		51	DQ5
144	As	544	DR10	DR11	52	DQ5 DQ7
145	AA	671	DR12	DR14	52	DQ5 DQ5
146	As	830	DR11	DR15	51,52	DQ6 DQ7
147	As	022	DR9	DR15	51,53	DQ5 DQ9
149	AA	531	DR7	DR13	52,53	DQ2 DQ7
150	AA	415	DR4	DR13	52,53	DQ6 DQ7
151	U	013	DR103			DQ5
152	AA	316	DR18		52	DQ4
154	C	342	DR4	DR13	52,53	DQ2 DQ6
157	AA	126	DR17		52	DQ2
161	As	730	DR10	DR12	52	DQ5 DQ7
162	AA	397	DR13	DR15	51,52	DQ6 DQ6
163	H	117	DR9	DR14	52,53	DQ5 DQ9
164	AA	088	DR8	DR18	52	DQ4 DQ6
166	C	307	DR8	DR17	52	DQ2 DQ4
167	U	839	DR4	DR4	53	DQ8
107		CON1	NEGATIVE (-) CONTROL BEAD			
111		CON2	NEGATIVE (-) CONTROL BEAD			
112		CON3	NEGATIVE (-) CONTROL BEAD			
108		IgG	POSITIVE (+) CONTROL BEAD			



Number of Beads Possessing Antigen

30 Beads for DR		
AG	#	%PRA
DR1	4	13%
DR103	1	3%
DR4	5	17%
DR7	3	10%
DR8	3	10%
DR9	4	13%
DR10	2	7%
DR11	5	17%
DR12	4	13%
DR13	6	20%
DR14	2	7%
DR15	4	13%
DR16	1	3%
DR17	5	17%
DR18	3	10%
DR51	5	17%
DR52	21	70%
DR53	12	40%
DQ2	10	33%
DQ4	4	13%
DQ5	11	37%
DQ6	8	27%
DQ7	11	37%
DQ8	2	7%
DQ9	2	7%

12 Beads for DQ		
DQ Enriched Beads		
AG	#	%PRA
DQ2	2	17%
DQ4	3	25%
DQ5	2	17%
DQ6	2	17%
DQ7	2	17%
DQ8	2	17%
DQ9	3	25%

BEAD	DONOR	DRB1*		DRB3*		DRB4*		DRB5*		DQB1*		DQA1*	
106	814	11:02	13:03	01:01	02:02					03:01		05:05	
113	373									02:02		02:01	
114	126									02:01		05:01	
115	031									04:02		04:01	
117	013									05:01		01:01	
118	634									06:01	06:02	01:02	01:03
119	699									03:01		05:05	
121	839									03:02		03:01	
122	012									03:03		03:02	
124	109									03:03	04:02	03:02	04:01G
125	461									03:03	05:01	01:01	02:01
126	484	03:01	11:01	02:02						02:01	03:19	05:01	05:05
127	467	08:04	15:03						01:01	03:19	06:02	01:02	04:01
128	106									03:02	06:01	01:03	04:01
129	560									03:01	04:02	04:01	05:05
130	373	07:01				01:01	01:03			02:02		02:01	
132	876	01:02	13:02	03:01						05:01	06:09	01:01	01:02
133	490	03:02	09:01	01:01		01:01				02:02	04:02	03:03	04:01
135	590	03:01	07:01	01:01		01:01				02:01	02:02	02:01	05:01
136	812	09:01				01:01	01:03			02:02		03:03	
137	328	01:01	04:01			01:01				03:01	05:01	01:01	03:03
138	653	01:01	11:01	02:02						05:01	06:02	01:01	01:02
139	941	01:02								05:01		01:01	
141	458	04:05	12:02	03:01		01:03				03:01	03:02	03:03	06:01
142	780	03:01	12:01	01:01	02:02					02:01	03:01	05:01	05:05
143	077	16:02						02:02	02:03	05:02		01:02	
144	544	10:01	11:01	02:02						03:01	05:01	01:01G	05:01G
145	671	12:01	14:01	01:01	02:02					05:01	05:03	01:04	01:05
146	830	11:06	15:01	02:02				01:01		03:01	06:02	01:02	05:05
147	022	09:01	15:01			01:03		01:01		05:02	03:03	01:02	03:02
149	531	07:01	13:03	01:01		01:01				02:02	03:01	02:01	05:05
150	415	04:07	13:01	01:01		01:03				03:01	06:03	01:03	03:03
151	013	01:03								05:01		01:01	
152	316	03:02		01:01						04:02		04:01G	04:04
154	342	04:05	13:01	01:01		01:03				02:02	06:03	01:03	03:03
157	126	03:01		02:02						02:01		05:01	
161	730	10:01	12:01	02:02						03:01	05:01	01:01G	05:05
162	397	13:02	15:03	03:01				01:01		06:02	06:09	01:02	
163	117	09:01	14:01	02:02		01:03				03:03	05:03	01:05	03:02
164	088	03:02	08:06	01:01						04:02	06:02	01:02	04:01
166	307	03:01	08:01	01:01						02:01	04:02	04:01	05:01
167	839	04:01	04:04			01:03				03:02		03:01	



LIFECODES Class II IDv2

LM2Q



See Instructions

CON1 (21 - 294)
CON2 (25 - 598)
CON3 (19 - 678)

EXPECTED RANGES

NAME

DATE

ID#

TECH

CalcCON

PRA

BDT

LOT 072811-LM2QR



2012-07-25

EXP: 25 JUL 12

BEAD	CON1 BAF	CON2 BAF	CON3 BAF	CalcCON BAF	SAMPLE	Sample/CON1	Sample/CON2	Sample/CON3	Sample/CalcCON	Sample/CON1-BAF	Sample/CON3 - BAF	Sample/CalcCON - BAF	INTRP
106	4.094	3.301	3.026	5.308									
113	4.742	3.444	3.332	6.306									
114	5.443	4.137	4.260	7.553									
115	4.263	3.177	3.431	6.667									
117	3.863	2.910	2.877	5.090									
118	4.127	3.523	3.236	5.570									
119	3.712	2.944	3.142	6.098									
121	4.218	3.144	3.118	6.067									
122	4.260	3.325	3.243	5.787									
124	4.483	3.575	3.575	6.796									
125	4.126	3.073	3.124	5.628									
126	3.818	3.263	2.881	5.190									
127	4.801	3.839	3.755	6.319									
128	4.563	3.514	3.334	5.959									
129	4.412	3.643	4.073	6.962									
130	3.563	2.706	2.914	5.059									
132	3.371	2.657	2.594	4.565									
133	3.674	2.999	2.761	5.080									
135	3.899	3.093	2.849	5.416									
136	3.808	3.103	2.996	5.154									
137	3.951	3.313	3.125	5.591									
138	3.309	2.625	2.647	4.541									
139	4.386	3.004	3.051	5.569									
141	3.900	3.095	2.857	5.100									
142	4.273	3.266	3.249	5.489									
143	4.811	3.482	3.602	6.158									
144	4.030	3.397	3.157	5.351									
145	4.206	3.408	3.094	5.444									
146	4.394	3.699	3.408	5.833									
147	3.946	3.206	3.082	5.401									
149	4.083	3.278	3.118	5.337									
150	4.088	3.382	3.047	5.323									
151	3.476	2.489	2.438	4.633									
152	3.733	2.851	2.695	5.025									
154	5.295	4.568	4.190	7.061									
157	3.287	2.436	2.485	4.412									
161	3.890	2.870	2.870	5.161									
162	5.440	4.428	3.941	7.456									
163	4.286	3.381	3.207	5.706									
164	4.897	3.881	3.747	6.642									
166	4.129	3.005	2.969	5.492									
167	5.543	4.214	4.059	7.508									