




DIACLONE Immunology Products

Anti-Human CD112 Azide Free

Clone: B-C12

CATALOGUE N° :	857.000.000
CLONE :	B-C12
ISOTYPE :	Mouse IgG2b
IMMUNISATION :	KG-1 cell line
HYBRIDOMA :	Myeloma X63/AG.8653 x Balb/c spleen cells
SPECIFICITY :	Recognises the Nectin-2, the polio retrovirus receptor-2 (PRR2)
WORKSHOP :	VI - VII
QUANTITY :	200 µg / 200 µl
PURIFICATION :	Ion exchange chromatography
FORMAT :	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
STORAGE :	Stable at +2-8°C for 12 months For longer storage, freeze aliquots
APPLICATION :	Suitable for flow cytometry Biological activity : not tested Immunoprecipitation : not tested Western Blot : not tested Suitable for immunohistochemistry
REFERENCES :	Lopez M. et al., J. Virol., 2000; 74(3) : 1267 - 1274 Lopez M. et al., Blood, 1998; 92 : 4602 - 4611



A typical staining pattern with the B-C12 monoclonal antibody of KG-1 cell line

Version 3 - 05.10

For research use only

20 YEARS OF EXPERTISE

With over 20 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

GLOBAL DISTRIBUTION

For the USA please contact our distributor Cell Sciences, Tel: 888-769-1246

email: info@cellsciences.com

For all other countries please visit our website for a full list of our distributors.

