




## DIACLONE Immunology Products

### Anti-Human CD13 Azide Free

Clone: B-F10

<b>CATALOGUE N° :</b>	854.100.000
<b>CLONE :</b>	B-F10
<b>ISOTYPE :</b>	Mouse IgG1
<b>IMMUNISATION :</b>	KG-1
<b>HYBRIDOMA :</b>	Myeloma X63/AG.8653 x Balb/c spleen cells
<b>SPECIFICITY :</b>	Recognises the aminopeptidase N, a 150-170 kDa protein
<b>QUANTITY :</b>	200 µg / 200 µl
<b>PURIFICATION :</b>	Ion exchange chromatography
<b>FORMAT :</b>	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
<b>STORAGE :</b>	Stable at +2-8°C for 12 months For longer storage, freeze aliquots
<b>APPLICATION :</b>	Suitable for flow cytometry Biological activity : not tested Immunoprecipitation : not tested Western blot : not tested Immunohistochemistry : not tested



A typical staining pattern with the B-F10 monoclonal antibody of monocytes

Version 3 - 05.10

For research use only

### 20 YEARS OF EXPERTISE

With over 20 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

#### GLOBAL DISTRIBUTION

For the USA please contact our distributor  
Cell Sciences, Tel: 888-769-1246

email: [info@cellsciences.com](mailto:info@cellsciences.com)

For all other countries please visit our website  
for a full list of our distributors.

