




DIACLONE Immunology Products


Anti-Human CD14 Unconjugated

Clone: B-A8

CATALOGUE N° :	854.113.020
CLONE :	B-A8
ISOTYPE :	Mouse IgG1
IMMUNISATION :	Monocytes
HYBRIDOMA :	Myeloma X63/AG.8653 x Balb/c spleen cells
SPECIFICITY :	Recognises the lipopolysaccharide receptor (LPS R), a 55 kDa protein
QUANTITY :	200 tests / 2 ml
PURIFICATION :	Ion exchange chromatography
FORMAT :	Phosphate-buffered saline with 1 % BSA and 0.1 % sodium azide
WORKING DILUTION :	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
STORAGE :	Stable at +2-8°C for 12 months For longer storage, freeze aliquots
APPLICATION :	Suitable for flow cytometry Immunoprecipitation : not tested Western blot : not tested Immunohistochemistry : not tested
REFERENCES :	Von Wolff M. et al., Mol. Hum. Reprod., 2002; 8(12) : 1096 - 1102



A typical staining pattern with the B-A8 monoclonal antibody of monocytes

 Sodium azide is a poisonous and hazardous substance, handle and dispose of properly

Version 1 - 05.01

For research use only

20 YEARS OF EXPERTISE

With over 20 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

GLOBAL DISTRIBUTION

For the USA please contact our distributor
Cell Sciences, Tel: 888-769-1246
email: info@cellsciences.com

For all other countries please visit our website
for a full list of our distributors.

