




DIACLONE Immunology Products

Anti-Human CD35 Unconjugated

Clone: UJ11

CATALOGUE N° :	853.403.020
CLONE :	UJ11
ISOTYPE :	Mouse IgG1,k
SPECIFICITY :	The CD35 antigen, CR1, is a complement component receptor with specificity for C3b, C4b and iC3b.
QUANTITY :	200µg / 200µl
PURIFICATION :	Affinity purification on protein G sepharose
FORMAT :	The monoclonal is presented in PBS with 15 mM sodium azide (pH 7.4)
WORKING DILUTION :	We recommend that each laboratory determine an optimum working titre for use in its particular application.
STORAGE :	Store at +2-8°C. For longer storage, freeze aliquots
APPLICATION :	Flow cytometry

 Sodium azide is a poisonous and hazardous substance, handle and dispose of properly

Version 2 – 01.07

For research use only

20 YEARS OF EXPERTISE

With over 20 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

GLOBAL DISTRIBUTION

For the USA please contact our distributor
Cell Sciences, Tel: 888-769-1246

email: info@cellsciences.com

For all other countries please visit our website
for a full list of our distributors.

