



## DIACLONE Immunology Products

### Anti-Human IL-12 p40 + p70 Azide Free

Clone: B-P24

|                       |   |
|-----------------------|---|
| <b>CATALOGUE N° :</b> | 855.130.005   |
| <b>CLONE :</b>        | B-P24   |
| <b>ISOTYPE :</b>      | Mouse IgG1  |
| <b>IMMUNISATION :</b> | Recombinant human IL-12   |
| <b>HYBRIDOMA :</b>    | Myeloma X63/AG.8653 x Balb/c spleen cells   |
| <b>SPECIFICITY :</b>  | Recognises both natural and recombinant p40 subunit alone and in its IL-12 form (associated with p35) and in its IL-23 form (associated with p19) |
| <b>QUANTITY :</b>     | 0.5 mg / 0.5 ml   |
| <b>PURIFICATION :</b> | Ion exchange chromatography   |
| <b>FORMAT :</b>       | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free  |
| <b>STORAGE :</b>      | Store at +2-8°C for 12 months<br>For longer storage, freeze aliquots  |
| <b>APPLICATION :</b>  | Blocks IL-12 induced proliferation on PHA and IL-2 activated lymphocytes<br>Suitable for intra-cytoplasmic flow cytometry                         |
| <b>REFERENCES :</b>   | Nagai, T. et al., J Immunol.,2003; 171(10): 5233-43.  |

Version 3 - 05.10

For research use only

### 20 YEARS OF EXPERTISE

With over 20 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacrone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacrone a clear choice to

### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

#### GLOBAL DISTRIBUTION

For the USA please contact our distributor

Cell Sciences, Tel: 888-769-1246

email: [info@cellsciences.com](mailto:info@cellsciences.com)

For all other countries please visit our website for a full list of our distributors.

