

p53 Antibody : B-P3 Azide free

ANTI HUMAN p53 AZIDE FREE

CATALOGUE N° : 853.040.005

CLONE : B-P3

ISOTYPE : Mouse IgG2a

IMMUNISATION : Recombinant p53 protein

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises wild type and mutant p53, cross-react with p53
in other primates
The epitope is within residues 18-30 of human p53

QUANTITY : 500 µg / 500 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm
and treated to remove endotoxins. Carrier and preservative
free

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Biological activity : not tested
Suitable for immunoprecipitation
Suitable for Western Blot
Suitable for immunohistochemistry in frozen and paraffin
sections
Suitable for intra-cytoplasmic flow cytometry

Version 1 – 05.01

For research use only